



## Integrated Triple Detection for the Analysis of Pectins by SEC

Greg Saunders  
Varian, Inc.

### Introduction

Pectins are a class of polysaccharide gum found naturally in fruits, such as apples, plums, grapes and cranberries. Structurally complex, pectins consist of 'smooth' and 'hairy' regions. The smooth regions are linear, partially methylated poly(D-galacturonic) acid, the hairy regions comprise alternating L-rhamnosyl and D-galacturonosyl residues containing L-arabinose and D-galactose branch points up to 20 residues long. As a result of this heterogenous nature, pectins adopt complex structures in solution. Applications of pectin are related to the formation of crosslinks through hydrogen bonding of the carboxylic acid groups, and include use as gelling agents, thickeners and water binders. Triple detection size exclusion chromatography employs a concentration detector, a viscometer and a light scattering detector to assess the molecular weight distribution and molecular structure of polymers without having to rely on column calibrations. This can be important when analyzing complex materials for which no structurally similar standards are available.

### Instrumentation

A sample of pectin was analyzed on the PL-GPC 50 Plus integrated GPC/SEC system from Varian, Inc. The instrument was operated at 30 °C and incorporated a refractive index detector, a PL-BV 400RT four capillary bridge viscometer and a PL-LS 15/90 light scattering detector (collecting scattered light at 15 ° and 90 °). Two PL aquagel-OH MIXED-H 8 μm columns were used for the analysis. These high performance columns offer excellent resolution over a very wide range of molecular weights, simplifying column selection and providing a versatile analytical system. The sample was prepared accurately in the eluent and filtered before injection through a 0.45 μm disposable filter. For the purpose of light scattering calculations, an average  $dn/dc$  value was used for the sample.

### Conditions

Sample: Pectin at 2 mg/mL  
Columns: 2 x PL aquagel-OH MIXED-H 8 μm, 300 x 7.5 mm  
Eluent: 0.2 M NaNO<sub>3</sub> + 0.01 M NaH<sub>2</sub>PO<sub>4</sub> adjusted to pH 7  
Flow Rate: 1.0 mL/min  
Injection Volume: 200 μL  
Sample Concentration: EasiVial™ PEO Standards, 0.1-0.5 mg/mL  
Polyacrylic Acid: approx 0.2 % w/v  
Detection: PL-GPC 50 Plus, RI, PL-BV 400RT, PL-LS 15/90

### Results and Discussion

Figure 1 shows an overlay of the triple detector chromatograms for the pectin sample. The chromatograms from the refractive index and light scattering detectors were clearly multimodal, as expected for a structurally heterogeneous material.




Figure 1. Triple detector chromatograms of pectin (autoscaled).




Figure 2. Molecular weight distribution plot calculated for the pectin.

From the viscometry and light scattering data, Mark Houwink (log intrinsic viscosity versus log M) and conformation (log radius of gyration versus log M) plots were generated for the pectin, shown overlaid in Figure 3.




Figure 3. Overlaid Mark Houwink and conformation plots for pectin.

The Mark-Houwink and to some extent the conformation plots show curvature over the entire molecular weight distribution, indicating a change in molecular density as a function of molecular weight, resulting from a variation in the relative amounts of 'smooth' and 'hairy' regions.

## Conclusion

The PL-GPC 50 Plus is a high resolution, cost effective integrated GPC system designed for operation from ambient to 50 °C. The standard system comprises precision solvent delivery, sample injection, high performance differential refractive index detection and a column oven, with fully integrated software control. When coupled with PL aquagel-OH MIXED-H columns, PL-BV 400RT viscometry and PL-LS 15 °/90 ° light scattering detectors, the PL-GPC 50 Plus makes maximum use of triple detection for the accurate determination of molecular weights of structurally complex and commercially important polymers.

*These data represent typical results. For further information, contact your local Varian Sales Office.*

*EasiVial, Varian and the Varian Logo are trademarks or registered trademarks of Varian, Inc. in the U.S. and other countries. © 2008 Varian, Inc.*

Varian, Inc.  
[www.varianinc.com](http://www.varianinc.com)  
 North America: 800.926.3000 – 925.939.2400  
 Europe: *The Netherlands*: 31.118.67.1000  
 Asia Pacific: *Australia*: 613.9560.7133  
 Latin America: *Brazil*: 55.11.3238.0400

