

Membrane Filter Unit

Membrane Filter Unit

**PRODUCT NAME:**

Nalgene® Analytical Filter Unit 130-4045

TEST TYPE:

Membrane filtration

PRODUCT DESCRIPTION:

This is an individually wrapped, pre-assembled, pre-sterilized, and disposable filter apparatus. The filter consists of a polypropylene upper chamber to pour the sample water into and a polystyrene lower chamber to collect the filtered water. The filter apparatus is meant to be for single use, but it can be re-used by sterilizing with an alcohol wipe (with an increased chance of cross-contamination). Filter apparatus includes upper and lower chambers with reference mark at 100 mL; removable, 47mm Triton™-free cellulose nitrate (CN) membranes, (0.45 µm pore size); cellulosic support pad; cellulosic vent plug side arm; quick-disconnect tubing adapter; and blue support plate.

MANUFACTURERS:

Nalge Nunc International Corporation
75 Panorama Creek Drive
Rochester, NY
14625, USA

web_master@nalgene.com

DISTRIBUTORS:

Various distributors in North America:

Cole-Parmer Canada
100-2500 Rue Jean-Perrin
Quebec, QC G2C 1X1 Canada
Toll Free: 800-363-5900

Phone: 514-355-6100

info@coleparmer.ca

www.coleparmer.ca

<https://www.coleparmer.com/i/thermo-scientific-nal-gene-130-4045-analytical-filter-units-0-45um-150-ml-12-pk/0673004>

ThermoFisher Scientific
3970 John's Creek Ct. Ste 500
Suwanee, GA
30024, USA

Tel: 1-800-955-5090 (In the United States)

Tel: +1-770-871-4725 (Outside the United States)

Fax: +1-770-871-4726

exports@thermofisher.com

AVAILABILITY:

USA and international

COST:

Contact ThermoFisher Scientific for pricing

OTHER NECESSARY PRODUCTS:

Vacuum pump



Nalgene® Analytical Filter Units (Source: Nalgene)

This document is open content. You are free to:

- Share – to copy, distribute and transmit this document
- Remix – to adapt this document

Under the following conditions:

- Attribution. You must give credit to CAWST as the original source of the document. Please include our website: www.cawst.org

CAWST and its directors, employees, contractors and volunteers do not assume any responsibility for and make no warranty with respect to the results that may be obtained from the use of the information provided.



B12, 6020 – 2 Street SE
Calgary, Alberta T2H 2L8
Canada
+1 (403) 243-3265
resources@cawst.org
cawst.org

To stay up-to-date and get support on:

- The latest updates to this document
- Other workshop and training related resources
- Using this document in your work

Visit resources.cawst.org