

Chemical Product Sheet

Water Sample Bag

**PRODUCT NAME:**

Whirl-Pak Bags

TEST TYPE:

Water sample bag

PRODUCT DESCRIPTION:

The Whirl-Pak bag is a sterile plastic bag that is sealed with a tab closure. The bags are meant for one-time use and are to be disposed afterwards. They are available in a variety of sizes ranging from 29 mL to over 5 L.

MANUFACTURER:

Nasco
901 Janesville Avenue
Fort Atkinson, WI
53538, USA

DISTRIBUTORS:

USA and Canadian Customers:
Tel: 1-800-558-9595
For all other international inquiries:
Tel: +1-920-563-2446
export@enasco.com
www.enasco.com/page/contact

AVAILABILITY:

US and international

COST:

Contact Nasco for pricing

STORAGE:


Should not be used for water temperatures above 82°C

OTHER NECESSARY PRODUCTS:

Permanent marker for writing label on bags




Whirl-Pak Bag (Source: Nasco)



Step 1

Step 2


Step 3



Step 4



Step 5



Step 6



Step 7

This document is open content. You are free to:

- Share – to copy, distribute and transmit this document
- Remix – to adapt this document



Under the following conditions:

- Attribution. You must give credit to CAWST as the original source of the document. Please include our website: www.cawst.org

CAWST and its directors, employees, contractors and volunteers do not assume any responsibility for and make no warranty with respect to the results that may be obtained from the use of the information provided.



B12, 6020 – 2 Street SE
 Calgary, Alberta T2H 2L8
 Canada
 +1 (403) 243-3265
resources@cawst.org
cawst.org

To stay up-to-date and get support on:

- The latest updates to this document
- Other workshop and training related resources
- Using this document in your work

Visit resources.cawst.org