

Chemical Product Sheet

Portable Incubator

**PRODUCT NAME:**

Portable Incubator (Product #2569900)

TEST TYPE:

Incubator for microbiological testing

PRODUCT DESCRIPTION:

This incubator can be used with a variety of power sources for reliable and accurate incubation. It has an operating range of 30 to 50°C, with an accuracy of +/-0.5°C at 37°C. The incubator can hold 6 Presence-Absence (P-A) test bottles, 39 most probable number (MPN) tubes, or 42 Petri dishes (50 mm each). Racks for P-A bottles, MPN tubes, and Petri dishes are sold separately. The external dimensions are 30.5 x 30.5 x 25.4 cm (12 x 12 x 10"). The power source is 12 volts direct current (vdc) using a vehicle cigarette lighter socket connection. Other power source options include 12 vdc rechargeable battery (12 hours of operation @ 37°C @ 20-25°C ambient) and 115/230 Vac Battery Eliminator. Includes power attachment for 12 vdc socket.

MANUFACTURERS:

Hach Company
Loveland, Colorado
USA

Tel: 1-800-227-4224

Fax: +1-970-669-2932

orders@hach.com

Outside the USA:

Tel: +1-970-669-3050

Fax: +1-970-461-3939

int@hach.com

www.hach.com

DISTRIBUTORS:

Various distributors worldwide

www.hach.com/global-distributor-support

AVAILABILITY:

USA and international

COST:

Contact Hach for pricing

MAINTENANCE:

Regular calibration, clean according to manufacturer's instructions

OTHER NECESSARY PRODUCTS:

Power source



Hach Portable Incubator
(Source: Hach)

This document is open content. You are free to:

- Share – to copy, distribute and transmit this document
- Remix – to adapt this document

Under the following conditions:

- Attribution. You must give credit to CAWST as the original source of the document. Please include our website:
www.cawst.org

CAWST and its directors, employees, contractors and volunteers do not assume any responsibility for and make no warranty with respect to the results that may be obtained from the use of the information provided.



B12, 6020 – 2 Street SE
Calgary, Alberta T2H 2L8
Canada
+1 (403) 243-3265
resources@cawst.org
cawst.org

To stay up-to-date and get support on:

- The latest updates to this document
- Other workshop and training related resources
- Using this document in your work

Visit resources.cawst.org